The iD-System™

The iD-Electrode™(non-sterile)

The iD-Paddle™(Sterile)

The iD-Smart Cable™(non-sterile)

User Manual

iD-Paddle™ “Large”

Ref.: iD-1200

Smart Clinical Products B.V

Ekkersrijt 4301

5692 DJ Son & Breugel

Netherlands

Tel: +31 499 465 770

Website: www.smartclinicalproducts.com

E-mail: info@smartclinicalproducts.com

© 2021 Smart Clinical Products B.V.

All rights reserved. Reproduction in whole or in part in any form or by any means, electrical, mechanical or otherwise is prohibited without the written consent of the copyright holder.

Manual version: SMART User Manual V10

Date of preparation: 10-Feb-2023

Symbols



Symbol	Description / Function	Reference
	Consult instructions for use	ISO 15223-1, 5.4.3
	Medical Device	ISO 15223-1, 5.7.7
	Indicates a single sterile barrier system	ISO 15223-1, 5.2.11
	Indicates a single sterile barrier system with protective packaging outside	ISO 15223-1, 5.2.14
	Caution: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1, 5.4.4
	Unique Device Identifier	ISO 15223-1, 5.7.10
	Do not re-use	ISO 15223-1, 5.4.2
	Do NOT re-sterilize	ISO 15223-1, 5.2.6
	Manufacturer information: This symbol is followed by the name and address of the device manufacturer.	ISO 15223-1, 5.1.1
	Indicates the date when the medical device was manufactured	ISO 15223-1, 5.1.3
	Indicates a medical device that has been sterilized using ethylene oxide	ISO 15223-1, 5.2.3
	This symbol means: Conforms to EC Directives and Regulations which apply to the device.	Article 20 Annex V
	Defibrillation-proof Type CF Applied part	IEC 60417-5336
	The WEEE symbol, indicating separate collection for WEEE- Waste of electrical and electronic equipment	DIRECTIVE 2012/19/EU, Annex IX
	Indicates the date after which the medical device is not to be used.	ISO 15223-1, 5.1.4
	Indicates the manufacturer’s catalogue number so that the medical device can be identified.	ISO 15223-1, 5.1.6
	Indicates that the medical device should not be used if the package has been damaged or opened	ISO 15223-1, 5.2.8
	Indicates the temperature limits to which the medical device can be safely stocked and exposed. Minimum +5°C or +41°F till Maximum +30°C or +86°F	ISO 15223-1, 5.3.7
	Indicates the temperature limits to which the medical device can be safely transported. Minimum -20°C or -4°F till Maximum +60°C or +140°F and this Max. 3days	
	Indicates the manufacturer’s batch code so that the batch or lot can be identified.	ISO 15223-1, 5.1.5

Table of Contents

1	General Description	4
1.1	Description of the Devices	4
1.2	Principle of Operation	5
1.3	Intended Purpose of the Devices	6
1.4	Clinical benefit	8
1.5	Claims for the Device.....	8
1.6	Performance Characteristics.....	8
1.7	Training requirements.....	8
1.8	Summary of Safety and Clinical Performance	8
1.9	Reporting of Incidents	8
2	Instructions For Use	9
2.1	Instructions For Use iD-System™	9
2.2	Instructions on Placement of the iD-Electrode™	10
2.3	Instructions on the use of the iD-Paddle™ and the iD-Smart Cable™.....	10

Chapter 1. General Description

1 General Description

1.1 Description of the Devices

The SMART Clinical Products’ iD-System™ comprises three separate devices:

1. The iD-Paddle™
2. The iD-Electrode™
3. The iD-Smart Cable™

The iD-Paddle is a disposable sterile medical device to be used as positive electrode for internal defibrillation during intra-thoracic surgery procedures. The iD-Paddle is surgically invasive, for transient use and intended to be used in contact with the patient’s heart to treat a cardiac condition; it is therefore classified as class III according to both the first and the third indent of rule 6.

The iD-Paddle comprises the ABS handle connected to the paddle’s head, also made of ABS and incorporating the electrode responsible for delivering the shock to the patient’s heart. Handle and head are joined via a flexible connection that allows the head to tilt in any direction facilitating the contact with the heart irrespective of the handle’s orientation. The “large” paddle has a contact area of 32cm²


The iD-Electrode is a disposable non-sterile device to be used as negative electrode for internal defibrillation in combination with the iD-Paddle. It is not invasive and is designed to be applied to the patient’s back (at the scapula level). The iD-Electrode is not an active device according to the MDR definition as it does not “act changing the density of or converting energy”. It is non-invasive and none of the conditions described in rule 2-4 are applicable to it. Therefore, it can be classified as class I according to rule 1.

The iD-Electrode comprises a 12,70cm x 6,98cm tin plate covered on one side with a conductive hydrogel (patient contacting side) and on the other side with a PE foam (non-patient contacting side). The tin plate is connected to an insulated wire and a connector that allows connection to the iD-Smart Cable.

The iD-Smart Cable is a non-sterile reusable accessory used to connect the defibrillator to the iD-Paddle and iD-Electrode. The iD-Smart Cable is, at this moment, compatible with different third-party defibrillators, as per the following list:

- iD-1500 – CU Medical’s LiFEGAIN CU-HD1, Progetti’s RESCUE Life, Bexen’s 700 and 800 REANIBEX
- iD-1510 – LIFEPAK 20, LIFEPAK 12 and LIFEPAK 15 of Physio-Control.
- iD-1530 – Philips DFM100 & Intrepid

The distal end of the iD-Cable, connecting with iD-Paddle and iD-Electrode, is the same for all models. The proximal end, connecting to the defibrillator, is customized to allow connection to the specific unit. The iD-Cable for connection with the Physio-Control’s defibrillator (iD-1510) also includes a shock button that controls the moment the energy is delivered to the patient’s heart.

The three parts of the iD-System: The iD-Paddle + The iD-Electrode + The iD-Smart Cable

Table 1 Device Description

Part # and description	Patient contact materials	Sterile	Contact with body / Location	Duration of use / Single- or multi-use
iD-1200 iD-Paddle “Large”	Aluminium (electrode) ABS (head/handle)	Y	Surgically invasive / Heart	Transient / Single use
iD-1100 iD-Electrode	Hydrogel KM 10G PE foam	N	Non-invasive / Intact skin	Transient / Single use
iD-1500-1510-1530 iD-Smart Cable	n/a	N	No patient contact	Transient / Multi-use

1.2 Principle of Operation

The iD-System operates by delivering a controlled electric shock to the patient’s heart. The energy is generated by a compatible defibrillator (not part of the system) and delivered via the iD-Paddle placed directly on the patient’s heart and the iD-Electrode applied to the patient’s back.

In the figure above, the position of the iD-Electrode on the back of the patient (Left), and the position of the iD-Paddle on the heart of the patient (right) are shown.

The iD-Paddle and iD-Electrode are connected to the defibrillator through the iD-Smart Cable.

The three components of the iD-System simply transmit the energy received from the defibrillator without transforming or converting it in any way.

Given the specific design of iD-Paddle and iD-Electrode, the procedure can be completed using one hand only, as opposed to the conventional paddles for internal defibrillation which require both hands to be operated.

1.3 Intended Purpose of the Devices

1.3.1. Intended Use

The iD-System is intended to be used together with a compatible defibrillator/monitor to deliver an electric shock to the surgically exposed heart.

1.3.2. Indications For Use

Internal defibrillation with the iD-System is indicated in patients experiencing ventricular fibrillation (VF), pulseless ventricular tachycardia (VT) or atrial fibrillation (AF) during open-chest surgery.

The iD-System is intended to be used by trained healthcare professionals in the OR (e.g. heart surgeons, anaesthesiologists, certified registered nurse anesthetists, etc.) or under supervision from such a person.

1.3.3. Target Population

The target population of the device under evaluation is the general population, with no specific exclusions.

1.3.4. Contra-indications

The iD-System is not indicated for external defibrillation.

iD-System is not indicated for patients <18 years old.

1.3.5. Warnings

General

- Do not exceed 50 joules while defibrillating
- Do not touch patient during defibrillation. Contact with the patient can lead to electric shock

iD-Paddle

- Do not use the iD-Paddle for the application of more than 20 shocks
- Do not use if package is open or damaged
- Do not use if the iD-Paddle shows any signs of damage
- Do not use without the iD-Electrode
- Do not reuse: risk of remnants of previous patient; other contamination; non-sterility; device might not function or not function correctly
- Do not re-sterilize

iD-Electrode

- Do not use if package is open or damaged
- Do not use if the iD-Electrode shows any signs of damage
- iD-Electrode is designed to be used together with the iD-Paddle only. Do not discharge any other hand-held paddle through this electrode
- Do not use additional gel on electrode
- Do not allow iD-Electrode to come into contact with any other materials such as ECG electrodes, wires, etc.
- Do not reuse: risk of contamination of the previous patient; device might not stick to the patient; device might not function correctly

1.3.6. Precautions

- Read instructions for use carefully before attempting to use.
- Do not open the iD-Paddle package until immediately prior to use
- Ensure the skin where the iD-Electrode is to be applied looks healthy, free from wounds or rashes
- Ensure the skin where the iD-Electrode is to be applied is dry, clean and free from cream or lotion. Dirt, grease or particles left on the skin may lodge themselves in the hydrogel of the iD-Electrode
- Do not use iD-Paddle or iD-Electrode beyond the expiry date shown on the labeling
- Ensure the iD-Electrode is properly applied to the skin. Inadequate patient contact may prevent the system to operate or may cause the current flow to concentrate with potential risk for skin burns
- Remove the iD-Electrode from the skin as soon as possible after use. Prolonged application to the patient’s skin can potentially lead to skin irritation
- Recommended storage temperature is +5°C to +30°C (+41°F to +86°F). Transport at temperatures from -20°C to +60°C (-4°F to +140°F) is limited to three days. When transport at these temperatures exceeds this time period, the shelf-life is reduced.

1.4 Clinical benefit

Ventricular fibrillation (VF), pulseless ventricular tachycardia (VT) and atrial fibrillation (AF) are potentially fatal occurrences during thoracic surgery. In those circumstances internal cardiac defibrillation and cardioversion are the only available options to restore a normal sinus rhythm and therefore increasing the chance of patient’s survival. The iD-System assists the connected defibrillator in achieving its intended purpose without having a direct therapeutic function itself.

1.5 Claims for the Device

The iD-System can transfer a shock pulse to the heart up to a maximum of 50 Joules.

The iD-Paddle’s 360° flexible head allows to easily achieve adequate contact with the heart.

The iD-System is designed to be an easy-to-use device operated with a single hand and to prevent misuse deriving from complex handling.

1.6 Performance Characteristics

Technical Characteristics	Related standard
Maximum 20 shots	-
Connection cable; specific per defibrillator mark and model	-
Cable length: ≥ 2 meter	IEC 60601-2-4
Large model surface area: ≥ 32 cm ²	IEC 60601-2-4
Small model surface area: ≥ 9 cm ²	IEC 60601-2-4
Shelf life ≥ 3 years	-
Pull strength on wire: ≥ 30 N	-
EO sterilization, sterile product	ISO 11135
Pouch packaging; pouch has no holes; seal strength > 1.2 N/mm ²	ISO 11607-1 ISO 11607-2
Biocompatible	ISO 10993-1 ISO 10993-5 ISO 10993-10 ISO 10993-7
Vibration and shock resistant	ASTM D4169
Storage temperature: +5°C to +30°C (+41°F to +86°F) Operating temperature: 0°C to +40°C (32°F to +104°F) Transport temperature: -20°C to +60°C (-4°F +140°F) and this Max. 3days	-
Meet applicable clauses of the electrical safety standard (product is CF applied part)	IEC 60601-1
Meet applicable specific defibrillator standard	IEC 60601-2-4

1.7 Training requirements

The intended users of the iD-System are medically trained healthcare professionals such as a Heart surgeon.

All users should be trained in the use of defibrillation of the heart.

All users are expected to have read this manual.

1.8 Summary of Safety and Clinical Performance

The Summary of Safety and Clinical Performance (SSCP) is available at the website of the Manufacturer: <https://www.smartclinicalproducts.com/>









1.9 Reporting of Incidents

In case of incidents that have occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Chapter 2. Instructions For Use

2 Instructions For Use

2.1 Instructions For Use iD-System

1 Apply iD-Electrode™ 	2 Press the “On” Button ← 	3 Connect The iD-Cable™ ← 	
4 Check Joule protocol regarding Internal defibrillation. 	5 Connect iD-Paddle™ ← 	6 Connect iD-Electrode™ ← 	
7 Press the “Charge” Button ← 	8 Press iD-Paddle on the heart 	9 Press the “shock” Button 	10 Repeat steps 7, 8 and 9 as needed. 

1: iD-Paddle


2: iD-Electrode

3: iD-Cable

2.2 Instructions on Placement of the iD-Electrode

1. Clean and dry skin site and shave if required.
2. Remove electrode from package.
3. Remove the electrode protective liner.
4. Position electrode as pictured below. Do not place on broken skin.
5. Smooth the electrode from the center outward to the edges with fingertips to ensure that there are no air pockets between the gel and the patient’s skin.
6. Connect the electrode to the defibrillator.
7. Follow the defibrillator instructions for use.
8. After use, remove electrode by slowly peeling pad from the patient’s skin.
9. Discard used electrode.

2.3 Instructions on the use of the iD-Paddle and the iD-Smart Cable.

No.	Instructions
1	Check expiry date before use and peel open the pouch
2	Press the “ON” button on the defibrillator.
3	Connect the iD-Smart Cable with the defibrillator (Each Brand defibrillator has his own iD-Smart Cable)
4	When the iD-Smart Cable is connected, confirm that the defibrillator’s energy level cannot go higher than 50 joules.
5	Connect the iD-Paddle with the iD-Smart Cable in the RED-Square hole.
6	Connect the iD-Electrode with The iD-Smart Cable in the BLACK-Round hole.
7	Select the energy level between 1 - 50 joules and press the Charge button.
8	Press the iD-Paddle down on the myocardium with sufficient pressure to ensure good contact to the heart
9	Press the shock button and a shock will be delivered on the myocardium of the heart.
10	Repeat steps 7,8 and 9 as needed.
11	At the end of operation, disconnect the iD-Paddle and discard.
12	The disposable devices described in this manual must be disposed of in compliance with the applicable local waste control regulations at the end of their service life. If you have questions regarding the disposal of the device, please contact the manufacturer or an authorized disposal company.

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



2.1 Instrucciones de uso del iD-System

1 Apply iD-Electrode™ 	2 Press the "On" Button ←  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	3 Connect The iD-Cable™ ←  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	
4 Check Joule protocol regarding Internal defibrillation.  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	5 Connect iD-Paddle™ ←  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	6 Connect iD-Electrode™ ←  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	
7 Press the "Charge" Button ←  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	8 Press iD-Paddle on the heart 	9 Press the "shock" Button  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	10 Repeat steps 7, 8 and 9 as needed. 

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



2 Mode d'emploi

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<p>1 Apply iD-Electrode™</p> 	<p>2 Press the "On" Button ←</p>  <p>Read Instructions Before Use</p>	<p>3 Connect The iD-Cable™ ←</p>  <p>Read Instructions Before Use</p>	
<p>4 Check Joule protocol regarding Internal defibrillation.</p>  <p>Read Instructions Before Use</p>	<p>5 Connect iD-Paddle™</p>  <p>Read Instructions Before Use</p>	<p>6 Connect iD-Electrode™</p>  <p>Read Instructions Before Use</p>	
<p>7 Press the "Charge" Button ←</p>  <p>Read Instructions Before Use</p>	<p>8 Press iD-Paddle on the heart</p> 	<p>9 Press the "shock" Button</p>  <p>Read Instructions Before Use</p>	<p>10</p> <p>Repeat steps 7, 8 and 9 as needed.</p> 

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



<p>1 Apply iD-Electrode™</p> 	<p>2 Press the "On" Button ←</p>  <p>Read Instructions Before Use</p>	<p>3 Connect The iD-Cable™ ←</p>  <p>Read Instructions Before Use</p>	
<p>4 Check Joule protocol regarding Internal defibrillation.</p>  <p>Read Instructions Before Use</p>	<p>5 Connect iD-Paddle™</p>  <p>Read Instructions Before Use</p>	<p>6 Connect iD-Electrode™</p>  <p>Read Instructions Before Use</p>	
<p>7 Press the "Charge" Button ←</p>  <p>Read Instructions Before Use</p>	<p>8 Press iD-Paddle on the heart</p> 	<p>9 Press the "shock" Button</p>  <p>Read Instructions Before Use</p>	<p>10 Repeat steps 7, 8 and 9 as needed.</p> 

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<p>1 Apply iD-Electrode™</p> 	<p>2 Press the "On" Button ←</p>  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	<p>3 Connect The iD-Cable™ ←</p>  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	
<p>4 Check Joule protocol regarding internal defibrillation.</p>  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	<p>5 Connect iD-Paddle™ ←</p>  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	<p>6 Connect iD-Electrode™</p>  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	
<p>7 Press the "Charge" Button ←</p>  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	<p>8 Press iD-Paddle on the heart</p> 	<p>9 Press the "shock" Button</p>  <p style="text-align: right; font-size: small;">⚠ Read Instructions Before Use</p>	<p>10 Repeat steps 7, 8 and 9 as needed.</p> 

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